


SC10[®] Benchtop Shutter Controller

SDK Manual



Copyright 2014-2023 Thorlabs, Inc.


Table of Contents

Chapter 1	Introduction	1
Chapter 2	C++ Software Development Kit	2
2.1.	sc10_cmd_library.h File Reference	2
	2.1.1. Functions	2
	2.1.2. Function Documentation	
Chapter 3	Python Software Development Kit	10
3.1.	SC10_COMMAND_LIB Namespace Reference	10
	3.1.1. SC10 Class Reference	
	3.1.2. Function Documentation	11
Chapter 4	LabVIEW Software Development Kit	15

SC10[®] SDK Manual Chapter 1: Introduction


Chapter 1 Introduction

User can start software development in C/C++ develop environment, LabVIEW and Python.



Chapter 2 C++ Software Development Kit

The user can start software development with Visual Studio 2017 or later versions. The supported files are in **\Thorlabs_SC10_C++SDK** under the **Sample** directory.



Copy **SC10CommandLib_win32.dll** or **SC10CommandLib_x64.dll** to your program folder, and make sure the library file(sc10_cmd_library.h) is in the same folder.

Below is the description of the header file sc10_cmd_library.h:

2.1. sc10_cmd_library.h File Reference

2.1.1. Functions

- **COMMANDLIB_API** int **List** (unsigned char *serialNo, int size) *list all the possible port on this computer.*
- **COMMANDLIB_API** int **Open** (char *serialNo, int nBaud, int timeout) open port function.
- COMMANDLIB_API int IsOpen (char *serialNo) check opened status of port
- **COMMANDLIB_API** int **Close** (int hdl) close current opend port
- **COMMANDLIB_API** int **Read** (int hdl, unsigned char *b, int limit)
- **COMMANDLIB API** int **Write** (int hdl, char *b, int size)
- **COMMANDLIB** API int Set (int hdl, char *c, int size)
- **COMMANDLIB_API** int **Get** (int hdl, char *c, unsigned char *d)
- **COMMANDLIB_API** int **Purge** (int hdl, int flag) *Purge the RX and TX buffer on port.*
- **COMMANDLIB_API** int **SetBaudRate** (int hdl, int value)
- **COMMANDLIB API** int **SetMode** (int hdl, int value)
- **COMMANDLIB_API** int **ToggleEnable** (int hdl)

- **COMMANDLIB_API** int **SetOpenTime** (int hdl, int value)
- **COMMANDLIB API** int **SetCloseTime** (int hdl, int value)
- **COMMANDLIB API** int **SetTriggerMode** (int hdl, int value)
- COMMANDLIB_API int SetExternalTriggerMode (int hdl, int value)
- **COMMANDLIB_API** int **SetRepeatCount** (int hdl, int value)
- **COMMANDLIB_API** int **GetBaudRate** (int hdl, int *value)
- **COMMANDLIB API** int **GetMode** (int hdl, int *value)
- **COMMANDLIB_API** int **GetEnableState** (int hdl, int *value)
- **COMMANDLIB_API** int **GetOpenTime** (int hdl, int *value)
- **COMMANDLIB API** int **GetCloseTime** (int hdl, int *value)
- **COMMANDLIB_API** int **GetTriggerMode** (int hdl, int *value)
- **COMMANDLIB_API** int **GetExternalTriggerMode** (int hdl, int *value)
- COMMANDLIB_API int GetRepeatCount (int hdl, int *value)
- COMMANDLIB API int GetClosedState (int hdl, int *value)
- **COMMANDLIB_API** int **GetInterlockTripped** (int hdl, int *value)
- **COMMANDLIB_API** int **GetStat** (int hdl, char *d)

Undocumented Command.

- COMMANDLIB_API int GetId (int hdl, char *d)
- COMMANDLIB_API int GetCommands (int hdl, char *d)
- **COMMANDLIB API** int **SaveSettings** (int hdl)
- **COMMANDLIB API** int **StoreConfiguration** (int hdl)
- **COMMANDLIB_API** int **LoadConfiguration** (int hdl)

2.1.2. Function Documentation

COMMANDLIB_API int Close (int hdl)

close current opend port

Parameters:

hdl	handle of port.

Returns:

0: success; negtive number : failed.

COMMANDLIB_API int GetBaudRate (int hdl, int * value)

Get baud rate

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	0: 9.6k, 1:115k

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetClosedState (int hdl, int * value)

Get closed state

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.

value	1: shutter is closed, 0: shutter is open

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetCloseTime (int hdl, int * value)

Get close duration

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	the shutter's close time in ms

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetEnableState (int hdl, int * value)

Get State

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
speed	0: the shutter is disabled, 1: enabled

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetExternalTriggerMode (int hdl, int * value)

Get Ex trigger mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	external trigger mode

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetId (int hdl, char * d)

get the SC10 id

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
d	output string (<255)

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetInterlockTripped (int hdl, int * value)

Get interlock tripped

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	1: interlock is tripped, otherwise 0

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetMode (int hdl, int * value)

Get operating mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	the mode value(1-5)

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetOpenTime (int hdl, int * value)

Get open duration

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	the shutter's open time in ms

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetPorts (char * serialNo)

list all the possible port on this computer.

Parameters:

serialNo	port list returned string include serial number and device descriptor, seperated
	by comma

Returns:

non-negtive number: number of device in the list; negtive number: failed.

COMMANDLIB_API int GetRepeatCount (int hdl, int * value)

Return the repeat count

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	repeat count, a value of 1-99

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int GetTriggerMode (int hdl, int * value)

Get trigger mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	0:internal trigger mode,1:external trigger mode

Returns:

0: success:

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int IsOpen (char * serialNo)

check opened status of port

Parameters:

serialNo	serial number of the device to be checked.

Returns:

0: port is not opened; 1: port is opened.

COMMANDLIB_API int List (char * serialNo)

list all the possible port on this computer.

Parameters:

serialNo	port list returned string include serial number and device descriptor, separated
	by comma

Returns:

non-negative number: number of device in the list; negative number: failed.

COMMANDLIB_API int LoadConfiguration (int hdl)

Load configuration from EEPROM

make sure the port was opened SUCCESSful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.

Returns:

0: SUCCESS; negative number: failed.

COMMANDLIB_API int Open (char * serialNo, int nBaud, int timeout)

open port function.

Parameters:

serialNo	serial number of the device to be opened, use GetPorts function to get exist list
	first.
nBaud	bit per second of port
timeout	set timeout value in (s)

Returns:

non-negtive number: hdl number returned successfully; negtive number: failed.

COMMANDLIB_API int SaveSettings (int hdl)

Save current baud rate and output trigger mode

make sure the port was opened SUCCESSful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

- 1	7 77	1 11 0
	hdl	handle of port.
	riui	number of port.

Returns:

0: SUCCESS; negative number: failed.

COMMANDLIB_API int SetBaudRate (int hdl, int value)

set SC10's serial baud rate

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	SC10 baud rate, 0 for 9.6k and 1 for 115k

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int SetCloseTime (int hdl, int value)

set close duration

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	shutter's close time in ms

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int SetExternalTriggerMode (int hdl, int value)

Set Ex Trigger mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	0:set the output trigger to follow the shutter output when the SH05 is
	connected, 1: force the trigger output to follow the controller output when an
	SH05 is equipped

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int SetMode (int hdl, int value)

Set operating mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	SC10 mode. 1-manual, 2-auto, 3-single, 4-repeat, 5-external gate

Returns:

0: success:

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int SetOpenTime (int hdl, int value)

set open duration

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.	
value	shutter's open time in ms	

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

COMMANDLIB_API int SetRepeatCount (int hdl, int value)

Set repeat count

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.	
value	set the repeat count when in repeat mode, a value of 1-99	

Returns:

0: success;

0xEA: CMD NOT DEFINED;

0xEB: time out;

COMMANDLIB_API int SetTimeout (int hdl, int timeout)

set SC10's timeout

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
timeout	timeout

Returns:

0: success;

COMMANDLIB_API int SetTriggerMode (int hdl, int value)

set the trigger mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.
value	0:internal trigger mode,1:external trigger mode

Returns:

0: success:

0xEA: CMD NOT DEFINED;

0xEB: time out;

COMMANDLIB_API int StoreConfiguration (int hdl)

Store configuration, save current settings(ex. mode, open time, close time) into EEPROM make sure the port was opened SUCCESSful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of port.

Returns:

0: SUCCESS; negative number: failed.

COMMANDLIB_API int ToggleEnable (int hdl)

Enable/Disable the shutter

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

Parameters:

hdl	handle of	port.
-----	-----------	-------

Returns:


0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

Chapter 3 Python Software Development Kit

The user can start software development with Python 3.7 or later versions. The supported files are in **\Thorlabs_SC10_PythonSDK** under the **Sample** directory.



Copy all the files to your python project, and make sure they are in the same folder. The SC10_COMMAND_LIB.py is the wrapper for SC10_COMMAND_LIB(SC10CommandLib_win32.dll / SC10CommandLib_x64.dll in C/C++ development environment). The SC10_COMMAND_LIB_EXAMPLE.py is the example code for how to use the python APIs.

User can import the SC10_COMMAND_LIB in 32-bit application with SC10CommandLib_win32.dll or 64-bit application with SC10CommandLib_x32.dll. If in 64-bit application, then you need to modify the __init__ fuction code " lib_path = "./ SC10CommandLib_win32.dll "" to " lib_path = "./ SC10CommandLib_x32.dll "" in SC10_COMMAND_LIB.py file.

3.1. SC10_COMMAND_LIB Namespace Reference

3.1.1. SC10 Class Reference

Public Member Functions

- def __init__ (self)
- def **open** (self, serialNo, nBaud, timeout)
- def **is_open** (self, serialNo)
- def close (self)
- def **set_baud_rate** (self, baud_rate)
- def **set_mode** (self, mode)
- def toggle_enable (self)
- def **set_open_time** (self, time)
- def **set_close_time** (self, time)
- def **set_trigger_mode** (self, mode)
- def set_external_trigger_mode (self, mode)
- def set_repeat_count (self, repeat_count)
- def **get_baud_rate** (self, baud_rate)

- def get mode (self, mode)
- def get_enable_state (self, enable_state)
- def **get_open_time** (self, time)
- def **get_close_time** (self, time)
- def get_trigger_mode (self, mode)
- def get external trigger mode (self, mode)
- def **get_repeat_count** (self, repeat_count)
- def **get_closed_state** (self, state)
- def get_interlock_tripped (self, tripped_state)
- def **get_id** (self, id)
- def save_settings (self)
- def store_configuration (self)
- def load_configuration (self)

Static Public Member Functions

- def list_devices ()
- def load_library (path)

Data Fields

hdl

Static Public Attributes

- None **sc10Lib** = None
- bool **isLoad** = False

3.1.2. Function Documentation

def close (self)

```
Close opened SC10 device
Returns:
0: Success; negative number: failed.
```

def get_baud_rate (self, baud_rate)

```
Get baud rate
Args:
baud_rate: 0: 9.6k, 1:115k
Returns:
0: Success; negative number: failed.
```

def get_close_time (self, time)

```
Get close duration
Args:
   time: the shutter's close time in ms
Returns:
   0: Success; negative number: failed.
```

def get_closed_state (self, state)

```
Get closed state
Args:
   state: 1: shutter is closed, 0: shutter is open
Returns:
   0: Success; negative number: failed.
```

def get_enable_state (self, enable_state)

```
Get State
Args:
    enable_state: 0: the shutter is disabled, 1: enabled
Returns:
    0: Success; negative number: failed.
```

def get_external_trigger_mode (self, mode)

```
Get Ex trigger mode
Args:
   mode: external trigger mode
Returns:
   0: Success; negative number: failed.
```

def get_id (self, id)

```
get the SC10 id
Args:
   id: output string (255)
Returns:
   0: Success; negative number: failed.
```

def get_interlock_tripped (self, tripped_state)

```
Get interlock tripped
Args:
    tripped_state: 1: interlock is tripped, otherwise 0
Returns:
    0: Success; negative number: failed.
```

def get_mode (self, mode)

```
Get operating mode
Args:
   mode: the mode value(1-5)
Returns:
   0: Success; negative number: failed.
```

def get_open_time (self, time)

```
Get open duration
Args:
   time: the shutter's open time in ms
Returns:
   0: Success; negative number: failed.
```

def get_repeat_count (self, repeat_count)

```
GReturn the repeat count
Args:
    repeat_count: repeat count, a value of 1-99
Returns:
    0: Success; negative number: failed.
```

def get_trigger_mode (self, mode)

```
Get trigger mode
Args:
   mode: 0:internal trigger mode,1:external trigger mode
Returns:
   0: Success; negative number: failed.
```

def is_open (self, serialNo)

```
Check opened status of SC10 device
Args:
    serialNo: serial number of SC10 device
Returns:
    0: SC10 device is not opened; 1: SC10 device is opened.
```

def list_devices ()[static]

```
List all connected mcm301 devices
Returns:
The mcm301 device list, each deice item is serialNumber/COM
```

def load_configuration (self)

```
Load configuration from EEPROM
Args:
Returns:
0: Success; negative number: failed.
```

def load_library (path)[static]

def open (self, serialNo, nBaud, timeout)

```
Open SC10 device
Args:
    serialNo: serial number of SC10 device
    nBaud: the bit per second of port
    timeout: set timeout value in (s)
Returns:
    non-negative number: hdl number returned Successful; negative number: failed.
```

def save_settings (self)

```
Save current baud rate and output trigger mode
Args:
Returns:
0: Success; negative number: failed.
```

def set_baud_rate (self, baud_rate)

```
set SC10's serial baud rate
Args:
   baud_rate: SC10 baud rate, 0 for 9.6k and 1 for 115k
Returns:
   0: Success; negative number: failed.
```

def set_close_time (self, time)

```
set close duration
```

```
Args:
   time: shutter's close time in ms
Returns:
   0: Success; negative number: failed.
```

def set_external_trigger_mode (self, mode)

def set_mode (self, mode)

```
Set operating mode
Args:
   mode: SC10 mode. 1-manual, 2-auto, 3-single, 4-repeat, 5-external gate
Returns:
   0: Success; negative number: failed.
```

def set_open_time (self, time)

```
set open duration
Args:
    time: shutter's open time in ms
Returns:
    0: Success; negative number: failed.
```

def set_repeat_count (self, repeat_count)

```
Set repeat count
Args:
    repeat_count: set the repeat count when in repeat mode, a value of 1-99
Returns:
    0: Success; negative number: failed.
```

def set_trigger_mode (self, mode)

```
set the trigger mode
Args:
   mode: 0:internal trigger mode,1:external trigger mode
Returns:
   0: Success; negative number: failed.
```

def store_configuration (self)

```
Store configuration, save current settings(ex. mode, open time, close time) into EEPROM Args:
Returns:
0: Success; negative number: failed.
```

def toggle_enable (self)


```
Enable/Disable the shutter
Args:
Returns:
0: Success; negative number: failed.
```

Chapter 4 LabVIEW Software Development Kit


The user can start software development with LabVIEW 2013 or later versions based on LabVIEW instrument driver mechanism. The supported files are in **\Thorlabs_SC10_LabVIEWSDK** under the **Sample** directory.

How to install


Copy to instr.lib folder under LabVIEW installation folder.



Destination folder: under %LabVIEW install path%\instr.lib Typically, C:\Program Files (x86)\National Instruments\LabVIEW 2013\instr.lib




How to find VI



How to use

1. From VI


Note: Before you open the SDK LabVIEW project, make sure the device has been connected to the computer.



2. From VI tree

Some classic data flow in VI tree.

Use the Example Finder to find examples demonstrating the usage of this instrument driver. To launch Example Finder, select "Find Examples..." from the LabVIEW Help menu.



3. From example

Examples show the classic usage. Examples' path: instr.lib\Thorlabs_SC10\Example

